# CERAMVEST<sup>®</sup> – Hi Speed

## Descripción

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## CERAMVEST<sup>®</sup> – Hi Speed

### Description

CERAMVEST-HI SPEED is an ultra fine, graphite-free, phosphate bound investment for use with all types of crown and bridges alloys.

### Working conditions:

As with all phosphatic coatings, the temperature of storage and preparation of the coating and of the mixing liquid is a determining factor in the expansion during setting, and therefore in the precision of fit and the cast surface.

The optimum preparation temperature is 20-22°C. Any variation in the optimum working temperature, whether of the powder, liquid or mixing utensils, or significant variations in the ambient temperature, affect the behavior of the product.

### Control of expansion:

The expansion can be regulated by diluting the mixing liquid EXPANSOR-B in different concentrations in order to obtain high precision with any type of alloy used. ( See Table 1.)

The maximum expansion is obtained by using the liquid EXPANSOR-B in a pure state (100%). To obtain lesser expansions, add distilled water

Note: Setting under pressure chamber reduces the setting expansion by approximately 10% (at 5 bar).

### Preparation of the liquid:

Important:  
CERAMVEST-HI SPEED must always be used with its special mixing liquid EXPANSOR-B, pure or diluted.

The dilutions can be effected by adding distilled or demineralised water, shaking slightly and then applying to the bottle an identifying label of the concentration made.

Concentrations of mixing liquid:

-1-

Type of alloy	Concentration of EXPANSOR-B
Precious	70%
Semi-precious	80-90%
Non-precious	80-100%

### Preparation of the cylinder:

CERAMVEST-HI SPEED can be used either with or without a metal ring. If a ring is used, place as usual a layer of insulation that will permit the expansion of the coating during setting. Then moisten slightly with water the insulation walls to prevent absorption of the coating.

In the case of plastic or silicone cylinders, apply a fine layer of Vaseline to facilitate subsequent extraction.

### Mixing proportions:

To obtain good results, it is essential to use the precise measurements.

For greater ease of use, CERAMVEST-HI SPEED is supplied in pre-dosed sachets of 150 gr., and a dosing tube is supplied for the liquid.

It is recommended always to use whole sachets and to measure the liquid with the original dosing tube.

In the case of using non-pre-dosed portions of powder, weigh the exact contents with a scale.

### The powder-liquid proportion is: 100 gr. Powder : 25 ml liquid

Size of cylinder	powder	mixing liquid
1x	60 gr.	15 ml
3x	1 sachet 150g.	37,5 ml
6x	2 sachets 150g.	75 ml
9x	3 sachets 150g.	112,5 ml

### Mixing:

Prepare the liquid EXPANSOR-B according to the concentration required (as indicated above), add the powder to the liquid and mix briskly by hand for 15-20 seconds to obtain a totally moist and uniform mass, then pass to mechanical vacuum mixing for 60 seconds depending on the mixer.

### Filling the cylinder:

CERAMVEST-HI SPEED flows over wax surfaces without the need to paint them with degreasing agents.

When the mixing time is completed, fill the cylinder, vibrating gently, for 10-15 seconds, then stop.

### Working time:

6-7 minutes ( 20-22 °C). The working time can be reduced when the liquid is at a higher temperature, with ambient temperatures above 22 °C, and with the use of contaminated utensils.

### Setting time:

Depending on the size of the cylinder, the setting time is 20-30 minutes

## Heating:

This investment allows two types of heating:

- **fast heating** : placing the cylinder in the furnace between the 20<sup>th</sup> and 30<sup>th</sup> minute after the mixing at the final casting temperature.
- **normal heating** : placing the cylinder in the furnace after the setting, making a drying process at 250°C for 40 minutes and following heat up to the casting temperature.  
This heating is recommended when the cylinder is introduced in the furnace after the 30<sup>th</sup> minute.

## Placing in the furnace:

- A few minutes before placing in the furnace scrape with a knife the upper surface of the cylinder.
- Place the cylinders in the furnace with the spouts downwards.
- Avoid contact of the cylinders with the base and walls of the furnace: use spacers or crown trays.
- In the case of several cylinders, also avoid contact between them.
- Due to the aggressive elimination of the wax, do not open the furnace during heating.

## Final temperatures:

700 - 750°C ( 1300-1380F ) for precious alloys
800 - 850°C ( 1475-1575F ) for semi-precious alloys
850- 900°C ( 1575-1650F ) for non-precious alloys.

Maintain the cylinder at the final temperature for:

Cylinder size	Heating time
1X	25 min
3X	35 min
6X	45 min

After the casting has been completed, bench-cool the ring to room temperature. Do not quench in cold water !!

## Storage and expiry

### CERAMVEST-HI SPEED powder

Store at 20-21 °C.

Keep the packages tightly closed in a dry place.

The coating is sensitive to moisture: prolonged exposure to the air may alter the characteristics of the product.

Expiry: 36 months after date of manufacture.

### EXPANSOR-B liquid

Store at 20-21 °C.

Prolonged exposure to temperatures lower than 5 °C or higher than 35 °C may produce freezing and/or destabilization of the product.

Keep the packages tightly closed after each use. Prolonged exposure to the air may cause crystallization.

Expiry: 24 months after date of manufacture.

Batch number and expiry date are marked on each package.

## Packaging

### CERAMVEST-HI SPEED

*Powder*

6669-420 4.2 kg ( 28 bags x 150 g )

*Powder + liquid kit*

6669-120 12 kg ( 8 x 1.5 kg ) + 3 x 1000ml

### EXPANSOR-B

6614-100 1000ml

6614-500 5000ml

## Physical characteristics: In accordance with ISO 9694 (1996)

Initial setting ( Vicat needle)	9 min.
Resistance to compression	4-6 MPa.
Setting expansion	2.0 %
Thermal expansion	1.2 %

Values for pure liquid EXPANSOR-B ( 100 % ) and at 21 °C

## Safety information

- The coating contains quartz powder. DO NOT INHALE THE POWDER! Continued inhalation of the powder may cause silicosis.
- During handling, avoid the formation of powder. Prevent inhalation with a protective mask, type P-2. Eliminate any spillages of powder with a damp medium.
- Avoid the formation of powder during demuffling by submerging the cylinder completely cold in water.
- Aspire and extract to the exterior the ammoniac vapors from the pre-heating.
- Avoid skin contact with the mixing liquid EXPANSOR-B. In the event of contact, wash with soap and water.
- In the event of contact with the eyes, wash under abundant running cold water and consult a doctor.

## Guarantee:

Our recommendations referring to the techniques applied, whether verbally, in writing or in the form of practical instructions, are based on our own experiences and/or tests, and must therefore be considered only as instructions.

For further information, see the safety data sheet of CERAMVEST-HI SPEED and EXPANSOR-B which can be requested from your usual supplier.

**protechno**

PRODUCTOS PROTÉSICOS S.L.

Polígono Empordà Internacional  
17469 VILAMALLA, GIRONA , SPAIN

Tel.: ( 34 ) 972 52 61 69

Fax: ( 34 ) 972 52 55 86

E-mail: [protechno@protechno.es](mailto:protechno@protechno.es)

[www.protechno.es](http://www.protechno.es)